
#### pational diagnostics



# NucliStain™

For Nucleic Acid Visualization

# INTRODUCTION

National Diagnostics offers an improved method of nucleic acid visualization with the development of NucliStain. NucliStain is a positive stain concentrate for the rapid detection of double and single stranded DNA and RNA in Agarose and Polyacrylamide gels.

NucliStain is intended as a replacement for the conventionally used ethidium bromide. Its sensitivity is nearly comparable to that of ethidium bromide, with the capability of detecting as little as 50 ng of DNA. However, the NucliStain advantages are that Nuclistain is much less hazardous than ethidium bromide and that the results are visible under normal lighting conditions therefore eliminating the need for expensive visualization equipment such as UV transilluminators, and for the time-consuming photographic recording of gels. As a result, there is no UV damage to the DNA.

NucliStain is extremely easy to use. Simply pour and dilute. The separated DNA appears as dark blue bands on a light blue background.

### METHOD OF USE

- 1. To prepare a working strength solution of NucliStain, add one part NucliStain concentrate to one-hundred parts of distilled water. Mix until the solution is uniform. Prepare sufficient stain solution to fully submerge the gel.
- 2. Immerse the gel in the stain solution for 20–30 minutes.
- 3. Immerse the gel in distilled water for 3–4 hours to destain the blue background. Bands can often be visualized minutes after destaining. However, for increased sensitivity, overnight destaining is recommended.

## ALTERNATE PROTOCOL

75-100 ng sensitivity, faster staining

- 1. Dilute 100  $\mu$ l NucliStain into 100 ml distilled water
- 2. Stain gel for 1-2 hours until bonds are visible. Light blue background can be reviewed by destaining in water.

#### ORDERING INFORMATION

NucliStain

EC-730

25 ml

#### For Additional Information And Order Placement:

United States	
Phone:	(800) 526-3867
Fax:	(404) 699-2077
Email:	info@nationaldiagnostics.com
	-
Europe	
Phone:	0115 982 1111
Fax:	0115 982 5275
Email:	slsinfo@scientificlabs.co.uk

Int'l © National Diagnostics Inc.

