

SAFETY DATA SHEET



Conforms to regulation (EC) no. EU 453/2010

SECTION 1 - IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1 Product Identifier

Product Name: AquaPor LM

Product Number: EC-204

1.2 Relevant Identified Uses of the Substance/Mixture and Uses Advised Against

Investigational research by professional users

1.3 Details of the Supplier of the Safety Data Sheet

Manufacturer

National Diagnostics
305 Patton Drive
Atlanta, GA 30036
(404) 699-2121
(800) 526-3867
info@nationaldiagnostics.com

1.4 Emergency Telephone Number

ChemTel Inc.

Contract number MIS8894340
1-800 255-3924 (United States, Canada, Puerto Rico & US Virgin Islands)
01-800-099-0731 (Mexico)
400-120-0751 (China)
000-800-100-4086 (India)
1-300-954-583 (Australia)
0-800-591-6042 (Brazil)
+1-813-255-3924 (All other regions)

SECTION 2 - HAZARDS IDENTIFICATION

2.1 Classification of the Substance or Mixture

Not a hazardous substance or mixture according to regulation (EC) No. 1272/2008.

2.2 Label Elements

This product has no labeling elements associated with EC directives or respective national laws.

2.3 Other Hazards

None found.

SECTION 3 - COMPOSITION/INFORMATION ON INGREDIENTS

3.2 Mixture

Chemical Names/Description

A polysaccharide complex.

Component List

Component	% Comp.	CAS #	EC #	1278/2008 Classification
Agarose	100	9012-36-6	232-731-8	N.A.

SECTION 4 - FIRST AID MEASURES

4.1 Description of First Aid Measures

Inhalation

Remove to fresh air. Get medical attention for any breathing difficulty.

Ingestion

Give several glasses of water to drink to dilute. If large amounts were swallowed, get medical advice.

Skin

Wash exposed area with soap and water. Get medical advice if irritation develops.

Eyes

Wash thoroughly with running water. Get medical advice if irritation develops.

4.2 Most Important Symptoms and Effects, Both Acute and Delayed

Inhalation

Coughing, sneezing.

Ingestion

Upset stomach, constipation.

Skin

No adverse effects expected.

Eyes

Itching, redness.

4.3 Indication of Any Immediate Medical Attention and Special Treatment Needed

Unknown/not applicable

SECTION 5 - FIRE FIGHTING MEASURES**5.1 Extinguishing media**

Water spray, dry chemical, alcohol-resistant foam, or carbon dioxide.

5.2 Special Hazards Arising from the Substance/Mixture**Hazardous Combustion Products**

Carbon monoxide and carbon dioxide.

Hazardous Decomposition Products

Carbon dioxide and carbon monoxide may form when heated to decomposition.

Hazardous Polymerization

Will not occur under normal conditions of use (See Sections 10.4 & 10.5).

5.3 Advice for Firefighters

In the event of a fire, wear full protective clothing and NIOSH-approved self-contained breathing apparatus with full facepiece operated in the pressure demand or other positive pressure mode.

5.4 Further Information

No data available.

SECTION 6 - ACCIDENTAL RELEASE MEASURES**6.1 Personal Precautions**

Wear appropriate protective equipment as specified in Section 8.

6.2 Environmental Precautions

Prevent discharge into the environment. Dike spills and stop leakage where practical. Do not allow material to enter drains.

6.3 Methods and Materials for Containment and Cleaning Up

Remove all sources of ignition. Ventilate area of leak or spill. Vacuuming or wet sweeping may be used to avoid dust dispersal. Small amounts of residue may be flushed to sewer with plenty of water.

6.4 References to Other Sections

For disposal information, see Section 13. For protective clothing and equipment, see Section 8.

SECTION 7 - HANDLING AND STORAGE**7.1 Precautions for Safe Handling**

Avoid contact and inhalation. Do not get in eyes, on skin, on clothing. Wash thoroughly after handling.

7.2 Conditions for Safe Storage (including any incompatibles)

Keep in a tightly closed container, stored in a cooled, dry, ventilated area away from sources of heat or ignition. Protect from physical damage. Isolate from incompatible materials (section 10).

Incompatibles

Oxidizers.

7.3 Specific End Uses

Investigational research by professional users

SECTION 8 - EXPOSURE CONTROLS/PERSONAL PRECAUTIONS**8.1 Control Parameters**

ACGIH Threshold Limit Value (TLV): 10 mg/m³

OSHA Permissible Exposure Limit (PEL): 15 mg/m³

8.2 Exposure Controls**Engineering Controls**

A system of local and/or general exhaust is recommended to keep employee exposures below the Airborne Exposure Limits. Local exhaust ventilation is generally preferred because it can control the emissions of the contaminant at its source.

Respiratory Protection

If exposure limits are exceeded, wear a full-face respirator with organic vapor cartridge and high efficiency dust mist filter. Beyond fifty times exposure limits or when exposure levels are not known, wear a full-face piece positive pressure respirator.

Eye Protection

Safety glasses. Maintain eye wash fountain and quick-drench facilities in work area.

Skin Protection

Wear protective gloves and clean body covering clothing.

SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on Basic Physical & Chemical Properties

a. Appearance	White powder	b. Odor	None
c. Odor Threshold	N.A.	d. pH	No data found.
e. Melting/Freezing Point (°C)	No data found.	f. Boiling point (°C)	N.A.
g. Flash Point (°C)	N.A.	h. Evaporation Rate	No data found.
i. Flammability	N.A.	j. Upper/Lower Flammability or Explosive Limits	N.A.
k. Vapor Pressure	No data found.	l. Vapor Density (Air = 1)	No data found.
m. Relative Density	1.7	n. Water Solubility	Moderate (1 to 10%)
o. Partition Coefficient n-octanol/water	N.A.	p. Autoignition Temperature (°C)	N.A.
q. Decomposition Temperature (°C)	N.A.	r. Viscosity	No data available.
s. Explosive Properties	N.A.	t. Oxidizing Properties	Not an oxidizer

SECTION 10 - STABILITY AND REACTIVITY

10.1 Reactivity

Not reactive under normal conditions of use.

10.2 Chemical Stability

Stable under recommended conditions of use and storage.

10.3 Possibility of Hazardous Reactions

Will not occur under normal conditions of use (See Sections 10.4 & 10.5).

10.4 Conditions to Avoid

No information found.

10.5 Incompatible Materials

Oxidizers.

10.6 Hazardous Decomposition Products

Carbon dioxide and carbon monoxide may form when heated to decomposition.

SECTION 11 - TOXICOLOGICAL INFORMATION

Product LD50 Values

Oral Rat LD50 (mg/kg)

Dermal Rabbit LD50 (mg/kg)

Component Cancer List Status

	NTP Carcinogen		IARC Category
	Known	Anticipated	
Agarose	No	No	None

Potential Health Effects

Inhalation

Nuisance dust. May cause coughing or sneezing if inhaled in large amounts.

Ingestion

No adverse effects expected. Very large ingestion may cause gastrointestinal upset due to gelation and bulking with possible constipation.

Skin

No adverse effects expected.

Eyes

No adverse effects expected but dust may cause mechanical irritation.

Carcinogenicity

Not listed as a known or anticipated carcinogen by NTP or IARC.

Mutagenicity

No information found.

Reproductive Toxicity

No information found.

Teratogenic Effects

No information found.

Routes of Entry

No information found.

Target Organ Statement

No information found.

SECTION 12 - ECOLOGICAL INFORMATION

12.1 Toxicity

	Vertebrates	Invertebrates	Algae	Microorganisms
Aquatic Toxicity (ppm unless otherwise noted)	No data	No data	No data	No data
	Birds	Arthropods	Plants	Microorganisms
Terrestrial Environment Toxicity (ppm unless otherwise noted)	No data	No data	No data	No data

12.2 Persistence and Degradability

No data

12.3 Bioaccumulative Potential

No data

12.4 Mobility in Soil

No data

12.5 Results of PBT and vPvB Assessment

No data

12.6 Other Adverse Effects

None

SECTION 13 - DISPOSAL CONSIDERATIONS

13.1 Waste Treatment Methods

Offer surplus or non-recyclable product to licensed disposal company. Disposal is subject to user compliance with applicable law and product characteristics at time of disposal. Dispose of packaging as product.

SECTION 14 - TRANSPORT INFORMATION

	ADR/RID	IATA	IMO	DOT
14.1 UN Number	N.A.	N.A.	N.A.	N.A.
14.2 Shipping Name	Not regulated.	Not regulated.	Not regulated.	Not regulated.
14.3 Hazard Class	N.A.	N.A.	N.A.	N.A.
14.4 Packing Group	N.A.	N.A.	N.A.	N.A.
14.5 Environmental Hazards	N.A.	N.A.	N.A.	N.A.
14.6 Special Precautions	N.A.	N.A.	N.A.	N.A.

SECTION 15 - REGULATORY INFORMATION

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance/Mixture

United States

TSCA Regulatory Statement

All intentional ingredients are listed on the TSCA Inventory.

SARA 311/312 Hazard Categories

Component	Fire	Pressure	Reactivity	Acute	Chronic
Agarose	No	No	No	No	No

Europe

EEC Regulatory

All intentional ingredients are listed on the European EINECS Inventory.

SECTION 16 - OTHER INFORMATION

Revisional Updates

- 4/26/2019 - Updated Section 1.4
- 5/29/2015 - Updated Sections 2.1 and 3.2
- 10/29/2013 - Released Version 1.0

NFPA Codes

Health 0 Flammability 0 Reactivity 0

Dangers

Agarose
None

MANUFACTURER DISCLAIMER: The information given herein is offered in good faith as accurate, but without guarantee. Conditions of the use and suitability of the product for particular uses are beyond our control. All risks of use of the product are therefore assumed by the user. Nothing is intended as a recommendation for uses which infringe valid patents or as extending license under valid patents. Appropriate warnings and safe handling procedures should be provided to handlers and users.