

SAFETY DATA SHEET



Conforms to regulation (EC) no. EU 453/2010

SECTION 1 - IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1 Product Identifier

Product Name: Basic Fuchsin

Product Number: HS-518

1.2 Relevant Identified Uses of the Substance/Mixture and Uses Advised Against

Investigational research by professional users

1.3 Details of the Supplier of the Safety Data Sheet

Manufacturer

National Diagnostics
305 Patton Drive
Atlanta, GA 30036
(404) 699-2121
(800) 526-3867
info@nationaldiagnostics.com

1.4 Emergency Telephone Number

ChemTel Inc.

Contract number MIS8894340
1-800 255-3924 (United States, Canada, Puerto Rico & US Virgin Islands)
01-800-099-0731 (Mexico)
400-120-0751 (China)
000-800-100-4086 (India)
1-300-954-583 (Australia)
0-800-591-6042 (Brazil)
+1-813-255-3924 (All other regions)

SECTION 2 - HAZARDS IDENTIFICATION

2.1 Classification of the Substance or Mixture

Classification according to Regulation (EC) No. 1272/2008 [EU-GHS/CLP]

H351 - Carcinogenicity (Category 2)

2.2 Label Elements

GHS LABEL ELEMENTS AND CLASSIFICATION

GHS Label Elements



WARNING

H351 - Suspected of causing cancer.
P202 - Do not handle until all safety precautions have been read and understood.
P281 - Use personal protective equipment as required.
P308+P313 - IF exposed or concerned: Call a POISON CENTER or doctor/physician.

2.3 Other Hazards

None found.

SECTION 3 - COMPOSITION/INFORMATION ON INGREDIENTS

3.1 Substance

Chemical Names/Description

Rosaniline chloride, Magenta O; Color index: 42510

Chemical Formula

$C_{19}H_{17}N_3 \cdot HCl$

Component List

Component	% Comp.	CAS #	EC #
Basic Fuchsin	> 99	632-99-5	211-189-6

SECTION 4 - FIRST AID MEASURES

4.1 Description of First Aid Measures

Inhalation

Remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Call a physician.

Ingestion

Induce vomiting immediately as directed by medical personnel. Never give anything by mouth to an unconscious person. Call a physician.

Skin

Immediately flush skin with plenty of soap and water for at least 15 minutes while removing contaminated clothing and shoes. Get medical attention. Wash clothing before reuse. Thoroughly clean shoes before reuse.

Eyes

Immediately flush eyes with plenty of water for at least fifteen minutes, lifting lower and upper eyelids occasionally. Get medical attention immediately.

4.2 Most Important Symptoms and Effects, Both Acute and Delayed**Inhalation**

No information found, but compound should be handled as a potential health hazard.

Ingestion

No information found, but compound should be handled as a potential health hazard.

Skin

May cause irritation, redness and pain.

Eyes

May cause irritation, redness and pain.

4.3 Indication of Any Immediate Medical Attention and Special Treatment Needed

Unknown/not applicable

SECTION 5 - FIRE FIGHTING MEASURES**5.1 Extinguishing media**

Use any means suitable for extinguishing surrounding fire.

5.2 Special Hazards Arising from the Substance/Mixture**Hazardous Combustion Products**

Burning may produce carbon monoxide, carbon dioxide, nitrogen oxides and hydrogen chloride.

Hazardous Decomposition Products

Burning may produce carbon monoxide, carbon dioxide, nitrogen oxides and hydrogen chloride.

Hazardous Polymerization

Will not occur under normal conditions of use (See Sections 10.4 & 10.5).

5.3 Advice for Firefighters

In the event of a fire, wear full protective clothing and NIOSH-approved self-contained breathing apparatus with full facepiece operated in the pressure demand or other positive pressure mode.

5.4 Further Information

No data available.

SECTION 6 - ACCIDENTAL RELEASE MEASURES**6.1 Personal Precautions**

Wear appropriate protective equipment as specified in Section 8.

6.2 Environmental Precautions

Prevent discharge into the environment. Dike spills and stop leakage where practical. Do not allow material to enter drains.

6.3 Methods and Materials for Containment and Cleaning Up

Ventilate area. Keep unnecessary and unprotected people away from area of spill. Pick up and place in a suitable container for reclamation or disposal, using a method that does not generate dust.

6.4 References to Other Sections

For disposal information, see Section 13. For protective clothing and equipment, see Section 8.

SECTION 7 - HANDLING AND STORAGE**7.1 Precautions for Safe Handling**

Avoid contact and inhalation. Do not get in eyes, on skin, on clothing. Wash thoroughly after handling.

7.2 Conditions for Safe Storage (including any incompatibles)

Keep in a tightly closed container, stored in a cooled, dry, ventilated area. Protect from physical damage. Isolate from incompatible materials (section 10).

Incompatibles

Strong oxidizers, acids

7.3 Specific End Uses

Investigational research by professional users

SECTION 8 - EXPOSURE CONTROLS/PERSONAL PRECAUTIONS

8.1 Control Parameters

ACGIH Threshold Limit Value (TLV): none established

OSHA Permissible Exposure Limit (PEL): None established

8.2 Exposure Controls

Engineering Controls

A system of local and/or general exhaust is recommended to keep employee exposures low. Local exhaust ventilation is generally preferred because it can control the emissions of the contaminant at its source.

Respiratory Protection

For conditions of use where exposure to the substance is apparent, consult an industrial hygienist. For emergencies, or instances where the exposure levels are not known, use a full-facepiece positive-pressure, air-supplied respirator.

Eye Protection

Safety glasses.

Skin Protection

Wear protective gloves and clean body covering clothing.

SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on Basic Physical & Chemical Properties

a. Appearance	Metallic green lustrous crystals.	b. Odor	Faint odor
c. Odor Threshold	N.A.	d. pH	No data
e. Melting/Freezing Point (°C)	250	f. Boiling point (°C)	No data
g. Flash Point (°C)	N/A	h. Evaporation Rate	No data
i. Flammability	N.A.	j. Upper/Lower Flammability or Explosive Limits	N/A
k. Vapor Pressure	No data	l. Vapor Density (Air = 1)	No data
m. Relative Density	1.22	n. Water Solubility	soluble
o. Partition Coefficient n-octanol/water	N.A.	p. Autoignition Temperature (°C)	N/A
q. Decomposition Temperature (°C)	N.A.	r. Viscosity	N.A.
s. Explosive Properties	N.A.	t. Oxidizing Properties	Not an oxidizer

SECTION 10 - STABILITY AND REACTIVITY

10.1 Reactivity

Not reactive under normal conditions of use.

10.2 Chemical Stability

Stable under ordinary conditions of use and storage.

10.3 Possibility of Hazardous Reactions

Will not occur under normal conditions of use (See Sections 10.4 & 10.5).

10.4 Conditions to Avoid

Incompatibles.

10.5 Incompatible Materials

Strong oxidizers, acids

10.6 Hazardous Decomposition Products

Burning may produce carbon monoxide, carbon dioxide, nitrogen oxides and hydrogen chloride.

SECTION 11 - TOXICOLOGICAL INFORMATION

Product LD50 Values

Oral Rat LD50 (mg/kg)

No data

Dermal Rabbit LD50 (mg/kg)

No data

Component Cancer List Status

	NTP Carcinogen		IARC Category
	Known	Anticipated	
Basic Fuchsin	No	No	None

Potential Health Effects**Inhalation**

No information found, but compound should be handled as a potential health hazard.

Ingestion

No information found, but compound should be handled as a potential health hazard.

Skin

May cause irritation, redness and pain.

Eyes

May cause irritation, redness and pain.

Carcinogenicity

No information found.

Mutagenicity

No information found.

Reproductive Toxicity

No information found.

Teratogenic Effects

No information found.

Routes of Entry

No information found.

Target Organ Statement

No information found.

SECTION 12 - ECOLOGICAL INFORMATION**12.1 Toxicity**

	Vertebrates	Invertebrates	Algae	Microorganisms
Aquatic Toxicity (ppm unless otherwise noted)	No data	No data	No data	No data
	Birds	Arthropods	Plants	Microorganisms
Terrestrial Environment Toxicity (ppm unless otherwise noted)	No data	No data	No data	No data

12.2 Persistence and Degradability

No data

12.3 Bioaccumulative Potential

No data

12.4 Mobility in Soil

No data

12.5 Results of PBT and vPvB Assessment

No data

12.6 Other Adverse Effects

None

SECTION 13 - DISPOSAL CONSIDERATIONS

13.1 Waste Treatment Methods

Offer surplus or non-recyclable product to licensed disposal company. Disposal is subject to user compliance with applicable law and product characteristics at time of disposal. Dispose of packaging as product.

SECTION 14 - TRANSPORT INFORMATION

	ADR/RID	IATA	IMO	DOT
14.1 UN Number	N/A	N/A	N/A	N/A
14.2 Shipping Name	Not regulated	Not regulated	Not regulated	Not regulated
14.3 Hazard Class	N/A	N/A	N/A	N/A
14.4 Packing Group	N/A	N/A	N/A	N/A
14.5 Environmental Hazards	N.A.	N.A.	N.A.	N.A.
14.6 Special Precautions	N.A.	N.A.	N.A.	N.A.

SECTION 15 - REGULATORY INFORMATION

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance/Mixture

United States

TSCA Regulatory Statement

All intentional ingredients are listed on the TSCA Inventory.

SARA 311/312 Hazard Categories

Component	Fire	Pressure	Reactivity	Acute	Chronic
Basic Fuchsin	No	No	No	Yes	No

Europe

EEC Regulatory

All intentional ingredients are listed on the European EINECS Inventory.

SECTION 16 - OTHER INFORMATION

Revisional Updates

4/26/2019 - Updated Section 1.4

5/29/2015 - Updated Sections 2.1 and 3.1

11/14/2013 - Released Version 1.0

NFPA Codes

Health 1 Flammability 0 Reactivity 0

MANUFACTURER DISCLAIMER: The information given herein is offered in good faith as accurate, but without guarantee. Conditions of the use and suitability of the product for particular uses are beyond our control. All risks of use of the product are therefore assumed by the user. Nothing is intended as a recommendation for uses which infringe valid patents or as extending license under valid patents. Appropriate warnings and safe handling procedures should be provided to handlers and users.